




# EDGE Method Development and Troubleshooting

# EDGE OneTouch Methods

OneTouch methods are a good starting point;  
however, some method development may  
be needed to achieve best results!



# EDGE Method Development

Create Method	
Settings	Solvent Water >
Cycles	Top Volume (mL) 15.0
Parameters	Bottom Volume (mL) 0.0
Wash	Rinse Solvent Water >
	Volume (mL) 0.0
	Temperature (°C) 100
	Hold Time 0:00 >

Service Mode 26.6 1:29 AM

Method Development

- The solvent and volume should be determined by the customer for an application
- Temperature, hold time, cycles and purge time may need some adjustments to achieve best results

# Solvent


EDGE is compatible with any solvent

Solvent	Polarity	Boiling Point (°C)	Vapor Pressure (kPA)	Application
Water	Most Polar	100	2.4	Nutraceuticals
Methanol		64.6	16.9	Plastics
Acetonitrile		81.65	9.9	Pesticide Residues
Acetone		56.2	30	Environmental
Dichloromethane		39.8	58	Environmental
Toluene		110.6	3.8	Dioxins
Petroleum Ether		30-60	31	Total Fat
Hexane	Least Polar	69	17.6	Environmental

We recommend using the same solvent for the extraction and rinse


# Volume

We recommend using a total extraction volume of at least 30 mL



- Rinse: Use for absorbent matrices
- Top Add: Make sure solvent covers the entire sample
- Bottom Add: Creates dispersive effect

# Temperature



# Hold Time

Increase hold time by 1 min increments to improve recoveries



When using a hold time you may have to increase the temperature offsets for 100 C and 175 C in settings in order to reach higher temperatures.

# Hold Time

**Settings**

System

Localization

Date/Time

Run

Solvents

Bottle Setup

**Calibration**

Heater Temperature (100°C)	4.0
Heater Temperature (175°C)	12.0
Power Ramp	30.0
Temperature	-0.5
Pressure	1.2
Syringe Top	0.0

Service Mode 24.0°C 9:42 AM


Increase if not reaching temperature during hold

Decrease if exceeding temperature during hold



# Cycles




Add Cycles for samples that have highly concentrated analytes






- If the total extraction volume for all cycles exceeds 40 mL then each cycle will be collected in a separate vial.
- For each additional vial used then one less sample will be able to be loaded


# Cycles

## Create Method




Settings	Cycle	Solvent	Top	Bottom	Rinse Solvent	Rinse Volume	Temp	Hold
Cycles	1	Water	15	0	Water	0	100	0:00
Parameters	2	Water	15	0	Water	0	100	0:00
Wash								



Service Mode 26.6 1:30 AM 

# Purge Time



# Purge Time

The screenshot shows the 'Settings' menu in 'Service Mode'. The 'Collection' section is expanded, showing the following parameters:

Parameter	Value
Drain Pressure	8
Drain Temperature	0
Max Drain Time	2:00
Purge Pressure	22
Max Purge Time	2:00

Annotations on the right side of the screen provide guidance:

- An upward-pointing arrow next to the value '8' is labeled 'If false possible clog error', indicating that increasing the Drain Pressure can help resolve this issue.
- Two downward-pointing arrows next to the values '2:00' (for Max Drain Time and Max Purge Time) are labeled 'Decrease is draining or purging too long', indicating that decreasing these times can resolve this issue.

The bottom status bar shows 'Service Mode 25.5°C 13:26'.

- False Possible Clog Error: Increase Drain Pressure
- Draining or purging too long after successful collection: Decrease Drain or Purge Time


# Possible Clog: Cooling

The screenshot displays the 'Settings' menu in 'Service Mode'. The 'Run' option is selected in the left-hand navigation menu. The 'Cooling' section is expanded, showing three parameters: 'Safe Temperature' set to 50, 'Safe Pressure' set to 20, and 'Sample Removal' set to 150. The values 50 and 20 are circled in black. The top status bar shows 'Service Mode 25.9°C 13:26'.

Section	Parameter	Value
Pressure Stability	Time	0:30 >
	Range	3
Cooling	Safe Temperature	50
	Safe Pressure	20
	Sample Removal	150

For Actuator to open the Safe Temperature and Safe Pressure must be reached




# Wash



## Carryover

- Increase volume of washes
- Increase number of washes
- Use different polarity solvents

# Wash

**Create Method**   

Settings





Cycles

Parameters

Wash



Wash	Solvent	Volume (mL)
1	Water >	10.0
2	Water >	10.0
3	Water >	10.0
4	Water >	10.0
5	Water >	10.0

5-30 mL

   **Service Mode 26.4 1:30 AM** 

# Calibrations

- All calibrations are done at CEM prior to shipping
  - Autosampler: X, Y and Z axis
  - Temperature: Band heater and internal thermocouple
  - Pressure: small internal pump used to purge sample
  - Volume: Syringe
- Do not recalibrate anything unless directed to do so





# Tips upon receiving your EDGE

- Currently all systems are run prior to shipping
  - There may be some water drops in the solvent lines
  - There may be some residual water on the system collected in the blank runs
  - Always run a blank run with each solvent line to prime the lines
  - The system will take longer to pre-heat for the first run on a cool system

# Troubleshooting

- **Volume recovery is low**



- Inspect the solvent lines for any bubbles and tighten fittings into rotary valve if bubbles present
- Run a blank with Q-Cup and Q-Disc to verify volume recovery
- Note that some solvent will be loss due to evaporation

- **Extract did not Drain**

- Verify of O-rings
- Run the flush to make sure the pressure gets below 10 psi
- Run a hot blank (water to 150 °C) to try and clear clog
- Run the clear clog procedure via instructions from bulletin

- **Common places to clog**

- Chamber
- Dispense needle
- Cooling coil




# Clearing a Clog


- Follow the instructions of the bulletin to determine where the clog is located
  - External Syringe to purge from dispense tip
  - Internal pump test for vent line and rigid tubing
  - Bottom add solvent test



# Clearing a Clog





# Clearing a Clog






# Clearing a Clog



# EDGE Tips



- Waste must be connected to run
- Verify calibration for position 4 and waste
- Run a blank to prime lines
- Verify enough solvent
- Verify a collection vial is present