



ORACLE

Better Moisture/Fat Analysis

Original Scope

Design Goal:

Create a universal fat system that removes the bottlenecks and limitations of reference chemistries and rapid techniques. Design a rapid system that no longer requires any form of method development.

Long-term vision:

To become the standard reference technique for fat testing worldwide


What is the ORACLE?

- First and only Universal Fat Analyzer
 - Rapid NMR that requires NO fat method development
- Accurately analyzes moisture and fat in ANY unknown food sample



How the ORACLE works

1. Send a radio frequency signal that interacts with the H^+ protons on the sample's fat molecules
2. Sends the protons to an excited state
3. Protons relax to their natural state, releasing energy which is equated to fat content



ORACLE vs other NMR

The ORACLE utilizes a breakthrough NMR technology developed by CEM that overcomes the deficiencies of previous NMR technologies.

Two improvements over Trac technology


1. Isolates detection of proton signals on fat molecules from all other sample components
2. Eliminates partial decay signals of varying fat molecules

Validation of Technology

- ~30 CRM's analyzed on ORACLE
 - Samples extensively tested in collaborative studies (typically 10+ certified laboratories)
- CEM outsourced 100's of samples to Eurofins and Silliker
 - Submitted samples in “blind” and “non-blind” fashion to capture true sample variability

Certificate of Analysis (COA)

Statement of measurement



Poultry feed – Proximates and Elements
Reference Material LGC7173

Assessed Values using Statutory Methods¹

| Constituent | Number of laboratories | Assessed value (g/100 g) | Uncertainty ³ (g/100 g) | Commission Directive | k value ³ |
|------------------|------------------------|--------------------------|------------------------------------|----------------------|----------------------|
| Moisture | 9 | 12.3 | 0.3 | 73/47/EEC | 2.23 |
| Oil ⁵ | 9 | 4.1 | 0.7 | 98/64/EC | 2.14 |
| Ash | 7 | 6.4 | 0.6 | 71/250/EEC | 2.23 |

Assessed Values²

| Constituent | Number of laboratories | Assessed value (g/100 g) | Uncertainty ³ (g/100 g) | Weight ⁴ (g) | k value ³ |
|------------------|------------------------|--------------------------|------------------------------------|-------------------------|----------------------|
| Moisture | 13 | 12.4 | 0.3 | 3 | 2.14 |
| Nitrogen | 10 | 2.56 | 0.19 | 0.5 | 2.31 |
| Oil ⁵ | 12 | 4.1 | 0.7 | 2.5 | 2.16 |
| Ash | 15 | 6.5 | 0.6 | 3 | 2.20 |
| Crude Fibre | 9 | 4.1 | 0.7 | 1 | 2.26 |
| Calcium | 10 | 1.44 | 0.15 | 1 | 2.11 |
| Chloride | 7 | 0.28 | 0.06 | 1 | 2.23 |
| Magnesium | 6 | 0.16 | 0.02 | 1 | 2.16 |
| Phosphorus | 10 | 0.63 | 0.03 | 1 | 2.11 |
| Potassium | 8 | 0.74 | 0.06 | 1 | 2.26 |
| Sodium | 9 | 0.17 | 0.05 | 1 | 2.16 |


| Constituent | Number of laboratories | Assessed value (mg/kg) | Uncertainty ³ (mg/kg) | Weight ⁴ (g) | k value ³ |
|-------------|------------------------|------------------------|----------------------------------|-------------------------|----------------------|
| Iron | 8 | 145 | 31 | 1 | 2.14 |
| Manganese | 8 | 131 | 19 | 1 | 2.18 |
| Zinc | 9 | 91 | 11 | 1 | 2.12 |

Notes:

1. These values have been assigned using only data derived from laboratories reporting analysis according to EEC method of analysis for the official control of feedstuffs⁶, as incorporated into UK law in 'The Feeding Stuffs (Sampling and Analysis) Regulations 1999'.
2. These values have been assigned using data derived from a variety of methods.
3. The uncertainty quoted is the half-width of the expanded uncertainty interval calculated using a coverage factor (k), providing a level of confidence of approximately 95 %.
4. Weight of sample taken for homogeneity assessment. These are the same for the analytes assessed using the statutory methods.
5. Determined according to Procedure B of 98/64/EC.
6. Sample treated under heating with hydrochloric acid before solvent extraction of the oil.

Date of Issue: November 2008
Amended: December 2008

Signed: _____
Gill Holcombe (Mrs)
for the Government Chemist




Material number: LGC7173
Batch number: 003
Page 1 of 7

Queens Road, Teddington, Middlesex, TW11 0LY, UK
Tel: +44 (0)20 8943 7000 • Fax: +44 (0)20 8943 2767 • www.lgcgroup.com

- Accompanies every CRM sample
- Information varies slightly based on where it was sourced (e.g. Muva Kempten vs. NIST)
- Assessed values for determined components (e.g. fat/oil, moisture, protein, etc)
- Explanation of Statistics
- Suggested sample sizes
- Handling/Preparation instructions
- Shelf life

CRM Samples Run on ORACLE


| Sample | Source | % Fat | +/- |
|---------------------------|--------------|-------|------|
| Wheat Flour | LGC | 1.39 | 0.17 |
| Processed Meat | LGC | 11.57 | 0.44 |
| Poultry Feed | LGC | 4.10 | 0.70 |
| Sweet Digestive Biscuit | LGC | 21.17 | 0.45 |
| Powdered Infant Formula | NIST | 30.43 | 0.95 |
| Cream Powder | Muva Kempten | 43.39 | 0.15 |
| Whole Milk Powder | Muva Kempten | 26.14 | 0.11 |
| Salted Butter | Muva Kempten | 82.00 | 0.78 |
| Butter | Muva Kempten | 82.43 | 0.12 |
| Butter | Muva Kempten | 82.62 | 0.82 |
| Fresh Cheese (Lact. Red.) | Muva Kempten | 24.19 | 0.38 |
| Processed Cheese | Muva Kempten | 14.58 | 0.18 |
| Fresh Cheese | Muva Kempten | 3.84 | 0.26 |
| Heavy Cream | Muva Kempten | 38.65 | 0.15 |
| Heavy Cream | Eurofins | 30.23 | 0.05 |
| Heavy Cream | Eurofins | 37.09 | 0.03 |
| Heavy Cream | Eurofins | 44.98 | 0.09 |
| UHT Milk | Muva Kempten | 1.71 | 0.01 |
| Milk (Past. Homog.) | Eurofins | 0.89 | 0.01 |
| Milk (Past Homog.) | Eurofins | 1.84 | 0.01 |
| Milk (Past Homog.) | Eurofins | 3.19 | 0.02 |
| Yogurt | Muva Kempten | 1.87 | 0.05 |
| Yogurt | Muva Kempten | 3.80 | 0.05 |
| Boiled Sausage | Muva Kempten | 19.75 | 0.36 |
| Heavy Cream | Eurofins | 36.83 | 0.11 |
| Parmesan Cheese | Muva Kempten | 24.98 | 0.12 |
| Milk Chocolate | Muva Kempten | 39.98 | 0.42 |



ORACLE Actalia Study

Further Validation from Actalia

- Actalia is a COFRAC accredited lab in France
 - Validates equipment for the dairy industry
 - Highly respected by ISO and IDF
 - Seen as “experts” in dairy analysis



Actalia Study

- 2 major conclusions from ORACLE testing
 1. The ORACLE “..reproducibility is lower than [better than] the reproducibility of the reference method.”
 2. The accuracy of the ORACLE compared to reference chemistry showed the “...regression slope (0.999) and the intercept (0.009) are not significantly different, respectively from 1.00 and zero (P=5%).”

Actalia Importance

- Accredited, respected third party company tested and approved the ORACLE as:
 - Accurate
 - Repeatable
 - Easier than Reference Chemistry
- This data shows ORACLE can replace both:
 - Reference Chemistry (mojonnier, gerber, etc)
 - NIR/FT-IR (no calibrations, accuracy = more \$\$)

ORACLE and AOAC

Is the ORACLE AOAC Approved???

- Short Answer: Yes!!
- Long Answer:
 - AOAC does not approve UNITS, only technology/methods
 - AOAC 2008.06 and PVM 1:2004 are approved methods using “Rapid determination of moisture/solids and fat in meat/dairy products by MW and NMR analysis”
 - ORACLE uses MW and NMR, so it is still an approved method

Fat and Moisture in Meat and Processed Meat

- Ground beef
- Chicken
- Turkey
- Pork
- All beef hot dogs
- Ham
- Pork sausage
- Potted meat

Fat and Moisture in Dairy Products

- Milk
- Cream
- Ice cream mix
- Yogurt
- Sour Cream
- Cheese
 - Mozzarella
 - Swiss
 - Cheddar
- Cream Cheese

Using the ORACLE

Two Ways to Operate

Rapid- SMART 6

- Process control labs that need rapid moisture & fat results
- Results = < 5 minutes
- Dry samples in the SMART 6 for moisture results and then analyze fat in ORACLE





High Throughput- Oven

- Testing labs running 50+ samples per day
- Dry samples overnight in oven
- Condition 1 hour in CEM Precision Heater Block and then analyze fat in ORACLE



SMART 6 + ORACLE Procedure



**SMART 6 Moisture Analysis
2-4 minutes**


**Condition in QuikPrep
<45 seconds**

**ORACLE Fat Analysis
30 seconds**


Air Oven Testing Sequence



**Dry in Oven
Overnight**



**Condition in Heater Block
30-60 minutes**



**ORACLE Fat Analysis
30 seconds**

ORACLE R&D Lab Benefits

- Removes the variability of reference testing, especially for blended products
- Gives users the ability to “blindly” test new or modified samples with confidence
- No need for revalidation of methods after product reformulations
- Differences in sample matrices have no effect on ORACLE accuracy

Global Repeatability

- All ORACLEs are designed to produce the same NMR signal
- Ensures consistent results across suppliers and manufacturers worldwide





ORACLE vs Everything

ORACLE Alternatives

- Reference Chemistry
 - Soxhlet, Mojonnier, Gerber, etc
- NIR/FT-IR
 - FOSS, Bruker, Perten
- SMART Trac
 - Old CEM Technology

Reference Chemistry

USDA Paper on Reference Testing

- USDA study on Reference Chemistry performance by Contract Labs
 - In general poor results, with 30-50% of results outside 1 St Dev
 - Shows the need to “question” reference chemistry

Reference materials to evaluate measurement systems for the nutrient composition of foods: results from USDA’s National Food and Nutrient Analysis Program (NFNAP)

| Class | Nutrient | Total CRMs | Total labs | Total values | Count of 0 to 1 | Count of 1 to 2 | Count of 2 to 3 | Count of > 3 | Percent > 2 | Percent > 3 |
|------------|-----------|------------|------------|--------------|------------------|--------------------|--------------------|--------------|--------------|-------------|
| Proximates | Moisture | 11 | 7 | 118 | 82 | 22 | 9 | 5 | 11.9 | 4.2 |
| | Protein | 9 | 5 | 106 | 60 | 24 | 12 | 10 | 20.8 | 9.4 |
| | Ash | 11 | 5 | 107 | 55 | 26 | 11 | 15 | 24.3 | 14.0 |
| | Total Fat | 11 | 6 | 129 | 52 | 39 | 15 | 23 | 29.5 | 17.8 |

Reference Chemistry Woes

| Milk Powder | | GC FAMES | Base Hydrolysis | Sohxlet | Roese Gottlieb | Blight and Dyer |
|---------------------------|------|----------|-----------------|---------|----------------|-----------------|
| Total Fat (g/100g sample) | Mean | 27.39 | 25.25 | 28.42 | 24.94 | 24.41 |
| | SD | 2.22 | 0.41 | 0.2 | 0.44 | 0.52 |
| | %RSD | 8.11 | 1.65 | 0.7 | 1.76 | 2.13 |

Aued-Pimentel et al. *Quim. Nova*, 2010, 33, 76 – 84

- Using the wrong reference extraction can lead to the wrong %Fat result
- “Mixed” samples containing various sample matrices are extremely difficult to extract properly
 - Never a problem for the ORACLE

Actalia Study

- 2 major conclusions from ORACLE testing
 1. The ORACLE “..reproducibility is lower than [better than] the reproducibility of the reference method.”
 2. The accuracy of the ORACLE compared to reference chemistry showed the “...regression slope (0.999) and the intercept (0.009) are not significantly different, respectively from 1.00 and zero (P=5%).”

Negatives of Chemical Extractions

- High Cost
 - Chemicals, disposal, labor, consumables, and more
- Safety Issues
 - Uses various hazardous solvents and exposed hot surfaces such as air ovens and hot plates
- Time per Test
 - Modified methods can take 15-20 minutes, full methods can take up to 16 hours
- Difficulty of SOP
 - Multiple opportunities for human error leading to poor repeatability and reproducibility

NIR and FT-IR

Moisture, Fat, SNF, and Protein

ORACLE Rapid Moisture/ Solids and Fat Analyzer

Sprint Rapid Protein Analyzer



- Better Accuracy than NIR and FT-IR technology
- No calibration maintenance or cost
- Typically 1 method for many products

Dairy Production Needs

CEM Recommended Equipment

- ORACLE (M/F), Sprint (P)



Products Tested

Milk
Cream
Liquid Whey
Cheese
Whey Powder
Milk Powder
Retentate
Condensed Milk
UF Milk
Additives
Ice Cream
Yogurt

FOSS Recommended Equipment

- FT120 (liquid) FoodScan (solid) NIRS DS2500 (powder)



*FoodScan can analyze powders but accuracy will be worse



Using AOAC to compare data

- AOAC data is unbiased, performed by certified laboratories
 - The best representation of true system accuracy
- CEM has Dairy AOAC approval for many products
 - FOSS only has AOAC approval for milk
 - FOSS has ISO approval for cheese (no useable data)
- CEM and FOSS both have AOAC studies for Meat products
 - Good comparison of system expectations

| | CEM | | | NIR | | |
|-------------|-------------------------------|-----------------------|----------------|--------------------------------|--------------------|----------------|
| | Moisture Analysis | | | | | |
| Sample Type | Reference Value (AOAC 950.46) | ORACLE (AOAC 2008.06) | Difference (%) | Reference Value (AOAC 950.46) | NIR (AOAC 2007.04) | Difference (%) |
| Beef | 67.31 | 67.07 | 0.24 | 65.23 | 62.30 | 2.93 |
| Pork | 60.07 | 60.05 | 0.02 | 61.17 | 60.51 | 0.66 |
| Chicken | 74.99 | 74.69 | 0.30 | 73.75 | 73.48 | 0.27 |
| Turkey | 74.67 | 74.39 | 0.28 | 73.85 | 73.69 | 0.16 |
| Hot Dog | 54.03 | 53.86 | 0.17 | 63.29 | 62.17 | 1.12 |
| | Average Difference | | 0.20% | | | 1.03% |
| | Fat Analysis | | | | | |
| Sample Type | Reference Value (AOAC 960.39) | ORACLE (AOAC 2008.06) | Difference (%) | Reference Value (AOAC 960.39) | NIR (AOAC 2007.04) | Difference (%) |
| Beef | 26.56 | 26.55 | 0.01 | 29.30 | 29.99 | 0.69 |
| Pork | 22.30 | 22.30 | 0.00 | 22.25 | 21.99 | 0.26 |
| Chicken | 2.91 | 2.88 | 0.03 | 3.17 | 3.25 | 0.08 |
| Turkey | 1.00 | 1.03 | 0.03 | 1.48 | 1.89 | 0.41 |
| Hot Dog | 29.79 | 29.85 | 0.06 | 15.39 | 15.05 | 0.34 |
| | Average Difference | | 0.03% | | | 0.36% |
| | Protein Analysis | | | | | |
| Sample Type | Reference Value (AOAC 981.10) | Sprint (AOAC 2011.04) | Difference (%) | Reference Value (AOAC 981.10) | NIR (AOAC 2007.04) | Difference (%) |
| Beef | 18.26 | 18.06 | 0.20 | 17.74 | 18.92 | 1.18 |
| Pork | 16.89 | 17.26 | 0.37 | 17.16 | 16.71 | 0.45 |
| Chicken | 21.73 | 22.25 | 0.52 | 22.36 | 22.74 | 0.38 |
| Turkey | 18.17 | 18.03 | 0.15 | 24.47 | 24.86 | 0.39 |
| Hot Dog | 9.41 | 9.80 | 0.39 | 16.42 | 15.25 | 1.17 |
| | Average Difference | | 0.29% | | | 0.94% |

CEM Dairy AOAC Data

| Sample Type | | CEM M/S% | CEM F% | Lab M/S% | Lab F% | AOAC M/S% | AOAC F% |
|-------------|---------|----------|--------|----------|--------|-----------|---------|
| Milk | Average | 45.60 | 39.93 | 45.50 | 39.94 | 45.57 | 39.93 |
| | St Dev | 0.04 | 0.08 | 0.07 | 0.08 | 0.00 | 0.04 |
| Heavy Cream | Average | 45.60 | 39.93 | 45.50 | 39.94 | 45.57 | 39.93 |
| | St Dev | 0.04 | 0.08 | 0.07 | 0.08 | 0.00 | 0.04 |
| Mozzarella | Average | 46.03 | 24.36 | 46.12 | 24.38 | 46.15 | 24.32 |
| | St Dev | 0.15 | 0.11 | 0.01 | 0.05 | 0.07 | 0.11 |
| Swiss | Average | 39.98 | 27.93 | 39.82 | 27.99 | 39.96 | 27.98 |
| | St Dev | 0.07 | 0.15 | 0.10 | 0.12 | 0.17 | 0.16 |
| Cheddar | Average | 36.68 | 31.32 | 36.76 | 31.29 | 36.76 | 31.29 |
| | St Dev | 0.12 | 0.11 | 0.10 | 0.14 | 0.05 | 0.13 |




CEM Proven as Accurate

- AOAC data shows that CEM is:
 - 5x more accurate for moisture (0.20% vs 1.03%)
 - 11x more accurate for fat (0.03% vs 0.36%)
 - 3x more accurate for protein (0.29% vs 0.94%)
- Data from AOAC and customers proves CEM is also more accurate for Dairy products
 - Better Accuracy = Better Products = Higher Cost Savings

Further Validation from Actalia

- Actalia is a COFRAC accredited lab in France
 - Validates equipment for the Dairy industry
 - Highly involved in ISO and IDF
 - Seen as “experts” in dairy analysis



NIR calculations

- Validation
 - Cost to build the method calibration
 - 20-50 samples X cost/test for reference X # of components X # of calibrations
 - 20-50 x \$50 x 3 = \$3,000-7,500/Calibration
 - Or buy premade calibration from FOSS \$5,000-7,000
- Annual Maintenance
 - Cost to use the system and maintain accuracy
 - 5 samples x cost/test for reference x # of components x frequency/year x # of calibrations
 - 5 x \$50 x 3 x 12 = \$9,000/Calibration

Cost of Ownership

| | NIR+FT-IR/FT-NIR | CEM |
|----------------|------------------|-----------------|
| Annual Tests | 5,000 | 5,000 |
| | Calibration Cost | Consumable Cost |
| 1 CALIBRATION | \$9,000 | \$6,000 |
| 4 CALIBRATIONS | \$36,000 | \$6,000 |
| 8 CALIBRATIONS | \$72,000 | \$6,000 |

Both prices approximated based on regional/volume pricing differences

- NIR costs based on suggested maintenance of ANN calibrations
 - 3 components (Moisture, Fat, Protein)
 - No reformulations or recalibrations, only typical maintenance
- CEM costs based on consumables for ORACLE
 - 2 components (Moisture, Fat)
 - List price (can be decreased based on purchase quantity)



FOSS ANN meat calibration - the key to reduced calibration costs



The ANN calibration has a huge advantage compared to other calibration techniques: A very robust calibration can be developed, with no limit as to how many samples can be included in the calibration.

With one ANN calibration it is possible to cover many different products, where you traditionally need to develop and maintain several calibrations. This means reduced calibration development and maintenance costs, as less reference analyses are required.

The purpose of this paper is to show how the ANN-calibration is a superior calibration method and a more cost-effective method compared to PLS.

The content is structured as follows:

- FoodScan measurement principle
- Calibration methods: PLS versus ANN
- The FoodScan ANN calibration for raw meat & meat products
- Case study

calibration – it comes with the FoodScan. It would only be necessary to verify each calibration (for slope & intercept adjustment of the ANN). The total reference analyses costs for verification of the ANN calibration in the example, when four constituents are determined, would be **4,800 US\$** (20 x 4 x 4 x15 US\$).

Likewise, the total annual costs of maintaining the 4 ANN versions would be **14,400 US\$** (12 x 5 x 4 x 4 x15 US\$).

Based on the example, the difference in number of reference analyses needed, and the related costs for operation of PLS and ANN can be summarised as shown in table 3.

Table 3. Comparison of PLS and ANN – based on example (table 1 and 2)

| | PLS | ANN | PLS | ANN | Savings (US\$) |
|--------------------------------------|---------------------|-----|------------|--------|----------------|
| | No. of ref. samples | | Costs (\$) | | |
| <i><u>Building calibrations:</u></i> | | | | | |
| - Actual building | 1400 | 0 | 72,000 | 0 | 72,000 |
| - Validation | 280 | 80 | 14,400 | 4,800 | 9,600 |
| <i><u>Annual maintenance</u></i> | | | | | |
| (monthly verification) | 840 | 240 | 43,200 | 14,400 | 28,800 |

5. Conclusion

- Both when establishing as well as in the annual maintenance of the calibrations, there are huge savings when using ANN.
- In many cases what the user saves by going from PLS calibrations to ANN can pay for the FoodScan within 1-2 years or even a shorter period.
- Particularly where new PLS calibrations have to be built, huge savings are achieved by using ANN instead.
- The higher the number of PLS calibrations used, the higher the savings by switching to ANN.

Dedicated Analytical Solutions

P/N No. 1025682
Issue No. 1
September 2003

FOSS Electric A/S
69 Slangenrupgade
DK-3400 Hilleroed
Company Reg.No. 7339 9815

Tel: +45 7010 3370
Fax: +45 7010 3371
E-mail: info@foss-electric.dk
Web: www.foss.dk

Conclusion vs NIR

1. Better Accuracy

- Proven by AOAC studies
- Better Accuracy = Better Process Control = More \$\$\$

2. Less \$\$\$ to maintain

- Consumables are cheaper than Calibrations

3. Easier to Use

- Fewer methods, no calibrations, easier sample prep

4. More Universal

- Same systems used for any product
 - Meat or Dairy, Liquid or Powder, etc

SMART Trac

Vs SMART Trac II

| | ORACLE | SMART Trac II |
|-------------------------------|-----------------------|-----------------------|
| Direct Analysis | Yes | Yes |
| Average Test Time | 3-4 minutes | 3-4 minutes |
| Method Development | No | Yes |
| Use with R&D/New Formulations | Yes | No |
| Footprint (W x D x H) | 15.6 x 22 x 14 | 14 x 14 x 22 |
| Consumables | 1-2 pads, 1 Trac film | 1-2 pads, 1 Trac film |

Method Development

SMART Trac II requires Method development

1. Samples being analyzed must first be tested via reference chemistry
2. Then raw NMR signals for each sample must be analyzed and plotted against their reference result
3. Samples are then separated into different methods based on which signals are linear
4. All future tests are then based on the linear calibration for that method

Where the ORACLE helps

- Removes the variability of reference testing present in Trac II method development
- Gives users the ability to “blindly” test new or R&D samples
- No need for revalidation of methods after product reformulations
- Reduces the need for sorting through methods to pick the right one
 - May still be present due to variance in moisture parameters, though to a lesser extent

ORACLE Customers

Global ORACLE Users

- Nestle Foods (#1 food/beverage company)
 - JBS (#5)
 - Tyson Foods (#6)
 - ADM (#7)
 - Cargill Meat (#9)
 - Kraft Heinz (#10)
 - Unilever (#12)
 - General Mills (#19)
 - Fonterra (#29)
 - Conagra (#30)
 - Eurofins
 - ALS Testing Labs
- ...and many more

Many chose CEM after having NIR in the past and losing money to bad accuracy & calibrations


ORACLE Success Stories

- ~200 ORACLEs sold in <2 years
 - Including Eurofins, Silliker, and other testing labs
- Nestle - 11 plants in Mexico
 - Passed validation for ALL Nestle Standards
 - Nestle hosted webinar in July to all Nestle factories
- Sigma Alimentos - 12 plants in Mexico
 - Success at Monterrey R&D, meeting with directors of Quality and Purchasing soon
- Schreiber, Lactalis, Saputo, Bel
 - All Cheese plants using ORACLE, all have plants in Mexico

Best Industries

- Processed Dairy (cheese, yogurt, ice cream, etc)
 - Need a fast, accurate test
 - Moisture, Fat, and Protein are all important
 - Need 2 NIR/FT-IR systems (more \$ than 2 CEM units)
- Meat Processing
 - Raw Meat - use ProFat for M/F/P in 1 cheap system
 - Cooked Meat - 1 ORACLE method for all products

ORACLE Sales Industries



QUESTIONS?

