



ISO 16756:2024 | IDF 259:2024
A New International Standard
for Rapid Fat Analysis of Milk and Milk Products

Impact of Fat Testing

Milk Production

- 965.7 million tons of milk produced annually (FAO Market Review 2023)
 - $\frac{1}{3}$ is used for further dairy processing
- Average QA test rate is about one test per 10,000 Liters
 - Total of ~40,000,000 tests each year globally



Current Cost

- Butterfat market rate \$3.6/lb
- More accurate testing has the potential to provide significant savings to the bottom line



Options

Reference Methods

Long Methods, Primary Methods

Pros

- Reliable
- Accurate
- Accredited

Cons

- Hazardous
- Difficult SOP
- Expensive



Rapid Methods

NIR, FTIR, TD-NMR

Pros

- Rapid
- Easy SOP
- Cheap Cost/Test

Cons

- Unreliable
- High Maintenance
- Indirect Analysis



Universal Fat Determination



Reliable

- No method calibrations
- No calibration maintenance

Accurate

- Results equal to or better than extraction and hydrolysis methods

Green

- Does not require solvents

Uses Accredited Technology

- ISO 16756:2024 | IDF 259:2024

ISO Method

ISO 16756:2024 | IDF 259:2024

Milk and milk products – Guidance for the application of Carr-Purcell-Meiboom-Gill (CPMG) pulsed time-domain nuclear magnetic resonance (TD-NMR) spectroscopy for fat determination

ORACLE is the only commercially available system that utilizes CPMG TD-NMR, meeting all requirements of ISO 16756:2024 | IDF 259:2024.



Workflow



Record sample weights and pre-dry samples



Prepare ORACLE samples



Condition samples in heater block



Place sample in ORACLE and press "Run"

The Proof is in the Pudding... Analysis


- Whole Wheat
- Wheat Flour
- Sweet-Digestive Biscuit
- Milk Chocolate
- Baking Chocolate
- Peanut Butter
- Boiled Sausage
- Meat Homogenate
- Processed Meat
- Poultry Feed
- Equine Feed
- Swine Feed
- Lamb Feed
- Dairy-Beef Feed
- Pheasant and Turkey



- UHT Milk
- Milk Powder
- Whole Milk Powder
- Infant Formula
- Powdered Infant Formula
- Cream Powder
- Salted Butter
- Butter
- Yogurt
- Fresh Cheese (multiple fat levels)
- Processed Cheese (multiple fat levels)
- Fresh Cheese (Lact. Red.)
- Parmesan Cheese




The Proof is in the Pudding... Analysis



- 36 CRMs initially validated on the ORACLE
- CEM outsourced 1000s of samples to contract labs
 - Submitted samples in “blind” and “non-blind” fashion to capture true sample variability

Actalia Cecalait Dairy Validation

- Experts in Dairy analysis
 - Validate and publish technologies for the dairy industry
- Independent study published by Actalia and available on cem.com
- Robust review of a wide variety of dairy products



Food Samples Tested

Over 2,000 unique sample types have been tested successfully

Dairy

raw, cultured, egg, alternative (soy, oat, nut, etc.), powders

Meat

raw, processed, smoked/cured, alternative (vegan)

Snack Foods

cookies, crackers, nuts, chips

Prepared Foods

condiments, frozen dinner, cereals, sauces, peanut butter

Alternative Foods

Plant-based meat & dairy, cultivated meat, and more

Confectionary

chocolate, syrup, cake, caramel, cocoa

Fats/Oils

rendered fats, cooking oils, fish oil, butterfats


Feeds/Meals

feed ingredients, additives, nutraceuticals, corn germ

More than 1000 units installed globally since release

- Manufacturing, R&D, Contract/Certified Labs, Government Labs, Academic

Universal Calibration – Large Food Testing Lab Validation



- Over 100 sample types
- Validated across 4x ISO 17025 accredited labs
 - 0-84% Fat range
 - RSD <1.0%
 - Bias is <0.03%

Options for Drying



Dry samples in oven



Prepare samples in ORACLE station



Place samples in heater block*



Analyze fat



Analyze moisture/solids and dry with SMART 6



Prepare samples in ORACLE station



* ORACLE can be paired with an optional robot for automated high-throughput sample testing



CPMG-NMR Benefits

Precise and reliable

- Results are comparable to reference chemistry
- Multiple independent laboratories have validated

Simple workflow integration

- Rapid results
- No harmful chemicals or waste
- Easy to use



Meets ISO 16756:2024 | IDF 259:2024 requirements for milk and milk products

